
Product Datasheet

Texas Red Conjugated Affinity Purified anti-Sheep IgG [H&L] [Donkey] (orb1675086)

Description

Texas Red Conjugated Affinity

[Donkey]

Species/Host

Donkey

Reactivity

Sheep

Conjugation

Texas Red

Tested Applications

FACS, FC

Immunogen

Sheep IgG whole molecule

Preservatives0.02 M Potassium Phosphate,
0.15 M Sodium Chloride, pH 7.2**Concentration**1.0 mg/ml (by UV absorbance
at 280 nm)**Storage**

For extended storage reconstitute product with 50% glycerol instead of water and then aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. Expiration date is one (1) year from date of restoration.

Note

For research use only

Application notes

Not specified. Optimal titers for applications should be determined by the researcher.

Expiration Date

12 months from date of receipt.