

Product Datasheet

Fluorescein Conjugated Affinity Purified anti-Donkey IgG [H&L] [Rabbit] (orb1675111)

Description	Fluorescein Conjugated Affinity Purified anti-Donkey IgG [H&L] [Rabbit]
Species/Host	Rabbit
Reactivity	Donkey
Conjugation	FITC
Tested Applications	FACS, FC
Immunogen	Donkey IgG whole molecule
Preservatives	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration	2.0 mg/ml (by UV absorbance at 280 nm)
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. Dilute only prior to immediate use. Expiration date is one (1) year from date of restoration.
Note	For research use only
Application notes	Not specified. Optimal titers for applications should be determined by the researcher.
Expiration Date	12 months from date of receipt.