

---

## Product Datasheet

### Texas Red Conjugated Affinity Purified anti-Sheep IgG [H&L] [Rabbit] (orb1675125)

**Description**

Texas Red Conjugated Affinity

[Rabbit]

**Species/Host**

Rabbit

**Reactivity**

Sheep

**Conjugation**

Texas Red

**Tested Applications**

FACS, FC

**Immunogen**

Sheep IgG whole molecule

**Preservatives**0.02 M Potassium Phosphate,  
0.15 M Sodium Chloride, pH 7.2**Concentration**1.0 mg/ml (by UV absorbance  
at 280 nm)**Storage**

For extended storage reconstitute product with 50% glycerol instead of water and then aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. Expiration date is one (1) year from date of restoration.

**Note**

For research use only

**Application notes**

Not specified. Optimal titers for applications should be determined by the researcher.

**Expiration Date**

12 months from date of receipt.