

## Product Datasheet

# Anti-PDGFR alpha phospho Y754 Antibody (orb1675140)

<b>Description</b>	Anti-PDGFR alpha phospho Y754 Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IHC, WB
<b>Immunogen</b>	Synthetic phospho peptide, containing Tyrosine 754 of Human PDGFR alpha (KEVSKYPSDIQRSL). The sequence is conserved in mouse and rat.
<b>Preservatives</b>	PBS, pH 7.4 with 0.05% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1ug/ul
<b>Storage</b>	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.