

Product Datasheet

Claudin 18.2 Antibody (orb1675186)

| | |
|----------------------------|--|
| Description | Claudin 18.2 Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC |
| Immunogen | Recombinant Protein |
| Form/Appearance | Tris Buffer, pH 7.3 - 7.7, with 1% rAlbumin and 0.1% Sodium Azide |
| Storage | This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt. |
| Note | For research use only |
| Application notes | 1:50-1:200 |
| Isotype | IgG |
| Clonality | Monoclonal |
| Expiration Date | 12 months from date of receipt. |