

Product Datasheet

Anti-MBP Tag Antibody (Chicken Polyclonal) (orb1675308)

Description	Anti-MBP Tag Antibody (Chicken Polyclonal)
Species/Host	Gallus
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgY
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.