
Product Datasheet

Anti-D Tag Antibody (Mouse Monoclonal) (orb1675311)

Description

Anti-D Tag Antibody (Mouse)

Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide (DYKDDDDK) coupled to KLH.
Form/Appearance	Liquid
Concentration	100ug/100uL
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use. Suggested starting dilutions are as follows: Western Blot at dilutions of 1:1000. Immunoprecipitation and immunostaining at dilutions of 1:1000.
Isotype	IgG2b with lambda light chain
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.