

## Product Datasheet

### Anti-CFP Tag Antibody (Monoclonal) (orb1675313)

<b>Description</b>	Anti-CFP Tag Antibody (Monoclonal)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Other
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IP, WB
<b>Immunogen</b>	GFP from the jellyfish <i>Aequorea victoria</i> N-Terminal peptide-KLH conjugates (first 20 aa)
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	100ug/100ul
<b>Storage</b>	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	WB 1:1000. For other applications, the researcher should determine the optimal working dilution.
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.