

Product Datasheet

Anti-Myc Tag Antibody (Mouse Monoclonal) (orb1675315)

Description	Anti-Myc Tag Antibody (Mouse Monoclonal)
Species/Host	Mouse
Conjugation	Unconjugated
Immunogen	A synthetic peptide corresponding to residues 410-419 (EQKLISEEDL) coupled to KLH.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	Western Blot (1:1000), Immunoprecipitation and Immunostaining (1:200). Other dilutions will need optimization by end user.
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.