


Product Datasheet

Human IL-13 Protein, GST Tag (orb1675324)

Description	Human IL-13 Protein, GST Tag is expressed from human 293 cells (HEK293). It contains AA Gly 21 - Asn...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-13
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-GST
Note	For research use only
Application notes	This protein carries a GST tag at the C-terminus. The protein has a calculated MW of 39 kDa. The protein migrates as 48-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Human IL-13 Protein, GST Tag (orb1675324) is expressed from human 293 cells (HEK293). It contains AA Gly 21 - Asn 132 (Accession # AAK53823.1).
MW	39 kDa
Uniprot ID	AAK53823.1
Expiration Date	6 months from date of receipt.



Immobilized
Human IL-13
Protein, GST
Tag...