
Product Datasheet

**Cynomolgus CA125 / MUC16 Protein, His Tag
(MALS verified) (orb1675329)**

Description

Cynomolgus CA125 Protein, His

cells (HEK293). It contains AA Pro 115...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

CA125

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 73.9 kDa. The protein migrates as 95-120 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

Pro 11544 - Pro 12188

Purity

90%

Source

Cynomolgus CA125 Protein, His Tag (orb1675329) is expressed from human 293 cells (HEK293). It contains AA Pro 11544 - Pro 12188 (Accession # XP_045234963.1).

MW

73.9 kDa

Expiration Date

6 months from date of receipt.