

---

## Product Datasheet

**Human Mesothelin / MSLN (391-598) Protein, Fc  
Tag (MALS verified) (orb1675337)**

**Description**

Human Mesothelin (391-598), Fc

293 cells (HEK293). It contains AA Ser 39...

**Reactivity**

Human

**Endotoxins**

0.2 EU per µg

**Conjugation**

Unconjugated

**Target**

Mesothelin

**Preservatives**

50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5

**Form/Appearance**

Powder

**Storage**

-20°C

**Tag**

N-Fc

**Note**

For research use only

**Application notes**

This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 49.6 kDa. The protein migrates as 50-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Protein Sequence**

Ser 391 - Ser 598

**Purity**

95%

**Source**

Human Mesothelin (391-598), Fc Tag (orb1675337) is expressed from human 293 cells (HEK293). It contains AA Ser 391 - Ser 598 (Accession # AAH09272.1).

**MW**

49.6 kDa

**Expiration Date**

6 months from date of receipt.