

---

## Product Datasheet

**Human GDNF / ATF / hGDNF Protein, premium grade (orb1675340)**

**Description**

Human GDNF Protein, premium grade is

contains AA Ser 78 ...

**Reactivity**

Human

**Endotoxins**

0.1 EU per µg

**Conjugation**

Unconjugated

**Target**

GDNF

**Preservatives**

PBS, pH7.4

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;  
-70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

**Tag**

Native

**Note**

For research use only

**Application notes**

This protein carries no "tag". The protein has a calculated MW of 15.1 kDa. The protein migrates as 15-22 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

90%

**Source**

Human GDNF Protein, premium grade (orb1675340) is expressed from human 293 cells (HEK293). It contains AA Ser 78 - Ile 211 (Accession # P39905-1).


**MW**

15.1 kDa

**Uniprot ID**
**P39905**
**Expiration Date**

6 months from date of receipt.

Human GDNF Protein, premium grade (250 µg of Human GDNF Protein, premium grade per well)


 Immobilized  
Human  
GDNF  
Protein,  
premium ...