


## Product Datasheet

### Anti-Human S100B (orb1675370)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Anti-Human S100B  |
| <b>Species/Host</b>        | Mouse   |
| <b>Reactivity</b>          | Human   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | WB  |
| <b>Immunogen</b>           | Recombinant human S100B protein (1-92/92 aa)  |
| <b>Target</b>              | S100B   |
| <b>Preservatives</b>       | 0.01M PBS(pH7.4)  |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG1  |
| <b>Clonality</b>           | Monoclonal  |
| <b>MW</b>                  | 10  |
| <b>Uniprot ID</b>          | <a href="#">P04271</a>  |
| <b>Dilution Range</b>      | WB=1:500-2000   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Western blot: Lane 1: Recombinant human S100B protein, N-His