

Product Datasheet

Mouse GPNMB/Osteoactivin ELISA Kit (orb1675402)

Description	Mouse GPNMB/Osteoactivin ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	<p>Detect Range: 156 pg/mL-10000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: GPNMB. Gene_full_name: Transmembrane glycoprotein NMB. Protein_name: Transmembrane glycoprotein NMB. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity macrophage inflamm. cancer tumor immunology tumor-associated antigens. Sequence_similarities: Belongs to the PMEL NMB family. Tissue_specificity: In early development mRNA is detected at highlevels in the outer layer of the retina. Later in developmentexpression gets restricted to the retinal pigment epithelium andiris. Overexpressed in mammary cancer cells possessing a strongbone metastatic phenotype.</p>
Expiration Date	Please enquire.