

Product Datasheet

Mouse Lymphocyte activation gene 3 protein ELISA Kit (orb1675564)

| | |
|----------------------------|--|
| Description | Mouse Lymphocyte activation gene 3 protein ELISA Kit |
| Reactivity | Mouse |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 0.15 ng/mL-10.0 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids. |
| Expiration Date | Please enquire. |