

Product Datasheet

Rat Luteinizing Hormone (LH) ELISA Kit (orb1675643)

Description Rat Luteinizing Hormone (LH) ELISA Kit

Reactivity Rat

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.468-30 ng/mL.
Sensitivity: 0.14 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.

Expiration Date Please enquire.