

## Product Datasheet

### Mouse IL-25(IL-17E) ELISA KIT (orb1675697)

<b>Description</b>	Mouse IL-25(IL-17E) ELISA KIT
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.03-2.0 ng/mL. Sensitivity: 15 pg/mL. Assay Time: 3 hour. Sample Volume: 20 µL. Sample Type: Cell culture supernatant, serum, plasma (EDTA, citrate, heparin). Detection method: Colorimetric.
<b>Expiration Date</b>	Please enquire.