

## Product Datasheet

# Plant hormone abscisic acid,ABA ELISA Kit (orb1675707)

<b>Description</b>	Plant hormone abscisic acid,ABA ELISA Kit
<b>Reactivity</b>	Plant
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.156 µg/mL-10 µg/mL. Sensitivity: 0.04 µg/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: plant tissues.
<b>Expiration Date</b>	Please enquire.