

Product Datasheet

Fish Estradiol(E2) ELISA kit (orb1675732)

Description	Fish Estradiol(E2) ELISA kit
Reactivity	Fish
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 40 pg/mL-1000 pg/mL. Sensitivity: 40 pg/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma, cell culture supernates, tissue homogenates.
Expiration Date	Please enquire.