

Product Datasheet

Bovine progesterone,PROG ELISA Kit (orb1675736)

Description	Bovine progesterone,PROG ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.15 ng/mL-70 ng/mL. Sensitivity: 0.2 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma.
Expiration Date	Please enquire.