

Product Datasheet

SARS-COV-2 Spike RBD(Omicron) Protein Antibody (orb1675788)

Description	SARS-COV-2 Spike RBD(Omicron) Protein Antibody
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of coronavirus SARS-COV-2 Spike RBD(Omicron) Protein (P0DTC2)(P0DTC2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35-40kDa
Uniprot ID	P0DTC2
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.