

Product Datasheet

BAK1 Rabbit Antibody (orb1675876)

Description	BAK1 Rabbit Antibody
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of arabidopsis thaliana BAK1(NP_567920.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68kDa
Uniprot ID	Q94F62
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.