

Product Datasheet

BAMBI Rabbit Antibody (orb1675936)

Description	BAMBI Rabbit Antibody
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human BAMBI(NP_036474.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55kDa
Uniprot ID	Q13145
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.