

Product Datasheet

GID1A Rabbit Antibody (orb1675994)

Description	GID1A Rabbit Antibody
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of arabidopsis thaliana GID1A(NP_187163.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38kDa
Uniprot ID	Q9MAA7
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.