

Product Datasheet

hnRNP A1 Rabbit Antibody (orb1676111)

Description	hnRNP A1 Rabbit Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human hnRNP A1(NP_002127.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29kDa/34kDa/38kDa
Uniprot ID	P09651
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.