

Product Datasheet

TRK fused gene (TFG) Rabbit Antibody (orb1676143)

| | |
|----------------------------|---|
| Description | TRK fused gene (TFG) Rabbit Antibody |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Recombinant fusion protein of human fused gene (TFG)(NP_006061.2). |
| Preservatives | PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 30kDa/31kDa/43kDa |
| Uniprot ID | Q92734 |
| Dilution Range | WB,1:500 - 1:2000 |
| Expiration Date | 12 months from date of receipt. |