

Product Datasheet

BES1 Rabbit Antibody (orb1676148)

| | |
|----------------------------|---|
| Description | BES1 Rabbit Antibody |
| Reactivity | A. thaliana |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | A synthetic peptide of arabidopsis thaliana BES1(NP_973865.1). |
| Preservatives | PBS with 0.05% proclin300,50% glycerol,pH7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 43KDa |
| Uniprot ID | Q9LN63 |
| Dilution Range | WB,1:500 - 1:2000 |
| Expiration Date | 12 months from date of receipt. |