



Product Datasheet

RPRD1B Antibody (orb1676254)


Description	RPRD1B Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived human RPRD1B recombinant protein (Position: R79-K272).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	Q9NQG5
Expiration Date	12 months from date of receipt.



IHC analysis of RPRD1B using anti-RPRD1B...



IHC analysis of RPRD1B using anti-RPRD1B...



IHC analysis of RPRD1B using anti-RPRD1B...