



Product Datasheet

Ogfr Antibody (orb1676300)

Description	Ogfr Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	E.coli-derived mouse Ogfr recombinant protein (Position: H217-K427).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Clonality	Polyclonal
MW	110 kDa
Uniprot ID	Q99PG2
Expiration Date	12 months from date of receipt.



Western blot analysis of Ogfr using anti-Ogfr antibody.



Flow Cytometry analysis of ANA-1 cells using anti-Ogfr antibody.