



Product Datasheet

RBMY1A1 Antibody (orb1676334)

Description	RBMY1A1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	E.coli-derived human RBMY1A1 recombinant protein (Position: R188-Y496).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51 kDa
Uniprot ID	P0DJD3
Expiration Date	12 months from date of receipt.



IF analysis of RBMY1A1 and Tubulin beta ...



Western blot analysis of RBMY1A1 using a...