


Product Datasheet

Human ICAM4 ELISA Kit (orb1676353)

Description	Human ICAM4 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A31-A240
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	Q14773
Sensitivity	15 pg/ml
Expiration Date	Please enquire.



Human
ICAM4
PicoKine
ELISA Kit
standard ...