

## Product Datasheet

# Woodchuck hepatitis B virus (isolate 8) (WHV) Protein X Antibody (orb1676457)

<b>Description</b>	Woodchuck hepatitis B virus (isolate 8) (WHV) Protein X Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	500 µg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P17401</a>
<b>Expiration Date</b>	12 months from date of receipt.