



Product Datasheet

Syntaxin 18/STX18 Antibody (orb1676524)


Description	Syntaxin 18/STX18 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	E.coli-derived human Syntaxin 18/STX18 recombinant protein (Position: M1-R313).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	Q9P2W9
Expiration Date	12 months from date of receipt.



IHC analysis of Syntaxin 18/STX18 using ...



IHC analysis of Syntaxin 18/STX18 using ...



IHC analysis of Syntaxin 18/STX18 using ...