


Product Datasheet

Mouse SDF-1/CXCL12 ELISA Kit (orb1676584)

Description	Mouse SDF-1/CXCL12 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: K22-M93
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin)
MW	13634 MW
Uniprot ID	H7BX38
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse SDF-1/CXCL12 PicoKine ELISA Kit st...