

Product Datasheet

Adenovirus Hexon Antibody (orb1676663)

Description	Adenovirus Hexon antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	CLIA, EIA, LF
Target	Adenovirus Hexon
Preservatives	0.05% Proclin 300
Form/Appearance	Supplied as a 0.22 µm filtered solution of PBS, pH 7.4
Storage	Short term (up to 7 days) store at 2–8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
Note	For research use only
Application notes	Method: Sandwich method. Use: Conjugate
Isotype	IgG
Clonality	Monoclonal
Purity	≥90% as determined by SDS-PAGE
Expiration Date	12 months from date of receipt.