

## Product Datasheet

### Adenovirus Hexon Antibody (orb1676664)

<b>Description</b>	Adenovirus Hexon antibody
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	CLIA, EIA, LF
<b>Target</b>	Adenovirus Hexon
<b>Preservatives</b>	0.05% Proclin 300
<b>Form/Appearance</b>	Supplied as a 0.22 µm filtered solution of PBS, pH 7.4
<b>Storage</b>	Short term (up to 7 days) store at 2–8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Method: Sandwich method. Use: Coating
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Expiration Date</b>	12 months from date of receipt.