

Product Datasheet

Recombinant Human TNF α Protein (orb1676739)

Description	Recombinant Human TNF α Protein
Tested Applications	CLIA, EIA, LF
Target	TNF α
Form/Appearance	Supplied as a 0.2 μ m filtered solution of PBS, pH 7.4
Concentration	Lot specific(\pm 20%)
Storage	Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
Note	For research use only
Application notes	Use: Calibration/quality control. Quality control: Mol Mass: 18.92 kDa AP Mol Mass:18.92 kDa Purity: \geq 90% as determined by reducing SDS-PAGE
Source	HEK293
Expiration Date	6 months from date of receipt.