

## Product Datasheet

# Recombinant Human Prolactin(PRL) Protein (orb1676751)

<b>Description</b>	Recombinant Human Prolactin(PRL) Protein
<b>Tested Applications</b>	CLIA, EIA, LF
<b>Target</b>	PRL
<b>Form/Appearance</b>	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4
<b>Concentration</b>	Lot specific(±20%)
<b>Storage</b>	Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Use: Calibration/quality control. Quality control: Mol Mass: 8.36 kDa AP Mol Mass: 8.50 kDa Purity: ≥90% as determined by reducing SDS-PAGE
<b>Source</b>	HEK293
<b>Expiration Date</b>	6 months from date of receipt.