


## Product Datasheet

### Cntn4 antibody (orb167682)

<b>Description</b>	Goat polyclonal antibody to Cntn4
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Canine, Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	Contactin 4 / Big-2 (mouse aa160-172)
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	WB: Approx 150kDa band observed in Mouse Brain lysates and not in the same lysates from the KO mice (calculated MW of 113kDa according to NP_001103219.1). Recommended concentration: 0.5-1.5µg/ml.
<b>Protein Sequence</b>	RRADGKPIARK
<b>Clonality</b>	Polyclonal
<b>MW</b>	113; 78.4
<b>Entrez</b>	<a href="#">152330</a>
<b>Dilution Range</b>	ELISA: 1:128000, WB: 0.5-1.5 µg/ml
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Mouse Olfactory...