

Product Datasheet

INDOL1 Polyclonal Antibody (orb1676848)

Description	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC, WB
Immunogen	KLH conjugated synthetic peptide derived from human INDOL1: 251-350/407
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 60 kDa
Uniprot ID	Q6ZQW0
Dilution Range	WB: 1:500-1:2000 ELISA: 1:5000-1:10000 IHC: 1:100-1:500 IF: 1:100-1:500 ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.