

## Product Datasheet

### ACACB antibody (orb167698)

<b>Description</b>	Goat polyclonal antibody to ACACB.
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Target</b>	ACACB (aa146-159)
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	WB: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 277kDa according to NP_001084.3. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring p
<b>Protein Sequence</b>	SPSKEDKKQANIKR
<b>Clonality</b>	Polyclonal
<b>MW</b>	277
<b>Entrez</b>	<a href="#">32</a>
<b>Dilution Range</b>	ELISA: 1:128000
<b>Expiration Date</b>	12 months from date of receipt.