

## Product Datasheet

### REG3 gamma monoclonal antibody (orb1677019)

<b>Description</b>	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	Recombinant fusion protein of human REG3 gamma. The exact sequence is proprietary.
<b>Preservatives</b>	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>MW</b>	~ 15, 19 kDa
<b>Uniprot ID</b>	<a href="#">Q6UW15</a>
<b>Dilution Range</b>	WB (1/500 - 1/1000), IHC (1/50 - 1/200)
<b>Expiration Date</b>	12 months from date of receipt.