

Product Datasheet

Alpha-internexin monoclonal antibody (orb1677042)

Description	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Alpha-internexin. The exact sequence is proprietary.
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 58 kDa
Uniprot ID	Q16352
Dilution Range	WB (1/500 - 1/1000), IHC (1/50 - 1/200), IF/ICC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.