

Product Datasheet

c-Met (Phospho-Y1234/Y1235) monoclonal antibody (orb1677313)

Description	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide of human c-Met (pY1234/Y1235). The exact sequence is proprietary.
Preservatives	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 140 kDa
Uniprot ID	P08581
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.