

Product Datasheet

TRK B (Phospho-Y705) Rabbit monoclonal antibody (orb1677735)

| | |
|----------------------------|--|
| Description | Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% BSA. |
| Species/Host | Rabbit |
| Reactivity | Human, Rat |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | A synthetic phosphopeptide corresponding to residues surrounding Tyr705 of human TrkB |
| Preservatives | Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin. |
| Concentration | 1mg/ml |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| MW | ~ 80,150 kDa |
| Uniprot ID | Q16620 |
| Dilution Range | WB (1/500 - 1/1000) |
| Expiration Date | 12 months from date of receipt. |