

Product Datasheet

Filamin A/B (Phospho-S2152/2107) Rabbit monoclonal antibody (orb1677738)

| | |
|----------------------------|--|
| Description | Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% BSA. |
| Species/Host | Rabbit |
| Reactivity | Human, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IF, IHC, WB |
| Immunogen | A synthetic phosphopeptide corresponding to residues surrounding Ser2152 of human Filamin A |
| Preservatives | Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin. |
| Concentration | 1mg/ml |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| MW | ~ 280 kDa |
| Uniprot ID | O75369 |
| Dilution Range | WB (1/500 - 1/1000), IHC (1/50 - 1/100), IF/ICC (1/50 - 1/100) |
| Expiration Date | 12 months from date of receipt. |