

Product Datasheet

SMAD3 (Phospho-S423/425) Rabbit monoclonal antibody (orb1677856)

Description	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% BSA.
Species/Host	Rabbit
Reactivity	Hamster, Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser423/Ser425 of human Smad3
Preservatives	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 52 kDa
Uniprot ID	P84022
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.