

## Product Datasheet

# COVID-19 S-clgG1 Cynomolgus monoclonal neutralizing antibody (orb1678008)

<b>Description</b>	1 mg/ml. Cynomolgus IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.4.
<b>Species/Host</b>	Monkey
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Preservatives</b>	1 mg/ml. Cynomolgus IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.4.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P0DTC2</a>
<b>Dilution Range</b>	ELISA: Use at an assay dependent dilution.
<b>Expiration Date</b>	12 months from date of receipt.