

Product Datasheet

COVID-19 S-mIgG1 Mouse monoclonal neutralizing antibody (orb1678009)

Description	1 mg/ml. Human IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.4.
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Preservatives	1 mg/ml. Human IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.4.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P0DTC2
Dilution Range	ELISA: Use at an assay dependent dilution.
Expiration Date	12 months from date of receipt.